

DID YOU KNOW?

A Pulse Oximeter may pose Potential Risks

*Risk Classification Under Medical Device Rules, 2017
Class B (Low to High Risk)*



Potential Risks

- Inaccurate Readings
- Reduced Accuracy with Poor Circulation
- External Interference such as dark skin pigmentation
- Skin Irritation

Precautions to be Taken

- Proper Device Use and Placement
- Monitor External Factors
- Double-Check with Clinical Symptoms
- Regular Maintenance and Cleaning

Any adverse event associated with these devices can be reported via following tools:-



HELPLINE
1800 180 3024

Available from: Monday to Friday
(9:00 AM to 5:30 PM)



MDAE
Reporting Form

Available on: www.ipc.gov.in



MDAE
Monitoring Centre

List Available on: www.ipc.gov.in



Adverse Drug
Reactions
Monitoring System

<https://adrmsipc.in/adrms/index.html>



Issued in Public Interest:

INDIAN PHARMACOPOEIA COMMISSION

National Coordination Centre, Materiovigilance Programme of India, Ministry of Health & Family Welfare, Govt. of India
Sector-23, Raj Nagar, Ghaziabad-201002

Email: lab.ipc@gov.in shatrunjay.ipc@gov.in Website: www.ipc.gov.in